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

Stage 1:	95℃	15min
Stage 2:	95℃	5s
	60℃	30s 30cycles
Stage 3:	Dissociation analysis	

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- [1] Specific and sensitive quantitative RT-PCR of miRNAs with DNA primers. BMC Biotechnol. 2011 Jun 25; 11: 70. PMC3135530.
- [2] Identification of circulating miRNA biomarkers based on global quantitative real-time PCR profiling. J Anim Sci Biotechnol. 2012 PMID:22958414.
- [3] Reliable reference miRNAs for quantitative gene expression analysis of stress responses in *Caenorhabditis elegans*. BMC Genomics. 2014. PMID:24656064.
- [4] Profiling microRNAs in lung tissue from pigs infected with *Actinobacillus pleuropneumoniae*. BMC Genomics. 2012. PMID:22953717.
- [5] The role of viral and host microRNAs in the Aujeszky's disease virus during the infection process. PLoS One. 2014. PMID:24475202.
- [6] MicroRNA polymorphisms and environmental smoke exposure as risk factors for oesophageal squamous cell carcinoma. PLoS One. 2013. PMID:24205249.
- [7] Plasma proANP and SDMA and microRNAs are associated with chronic mitral regurgitation in a pig model. Endocr Connect. 2013. PMID:24029364.
- [8] Advances in microRNA experimental approaches to study physiological regulation of gene products implicated in CNS disorders. Exp Neurol. 2012. PMID:22245616.
- [9] Maternally deposited germline piRNAs silence the transposon in somatic cells. EMBO Rep. 2013. PMID:23559065.
- [10] miRNA expression profile analysis in kidney of different porcine breeds. PLoS One. 2013. PMID:23372853
- [11] Determination of reference microRNAs for relative quantification in porcine tissues. PLoS One. 2012. PMID:22970213.