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

图 1, HaiGene TaqMan miRNAs 定量检测原理

## Material and Methods

The miRNA was extracted from the tissue or cell samples using Tissue&Cell microRNA Extraction Kit (HaiGene, Harbin, China) according to manufacturer's protocol. To investigate the miR-21 expression level in tissue samples, cDNA synthesis was carried out according to the protocol of TaqMan miRNA cDNA Synthesis Kit (HaiGene, Harbin, China). Briefly, the reaction mix consisting of 100 ng of miRNA, 1 µl of 5×TaqMan miRNA RT Solution A, was performed to add the polyA tail at 37 °C for 30min, 85 °C for 5min. Followed with adding 2 µl of 10×TaqMan miRNA RT Solution B, 2 µl of 10×TaqMan miRNA RT Primer, and 11 µl of ddH<sub>2</sub>O, which incubated at 30 °C 5min, 55 °C 60min followed by enzyme inactivation at 95 °C for 5 minutes. Quantification of miRNA was performed by HG TaqMan miRNA PCR Kit (HaiGene, Harbin, China). RealTime PCR was performed in 20 µl with 2 µl of cDNA product, 1 µl of miRNA TaqMan Assay, and 4 µl of 5×Golden HS TaqMan qPCR Mix under the following conditions: 95 °C for 15min, 40 cycles of 95 °C for 10s, and 60 °C for 30s in LightCycler 480 (Roche). MiR-21 levels were normalized to U6B RNA levels using the 2(-ΔCt) model. Real-time PCR assays were performed in duplicate for each sample, and the mean value was used for the calculation of miRNA expression levels. The statistical analyses were done using Microsoft Excel (Microsoft) both to calculate the SD and to test for statistically significant differences between the samples using a T test. A value P < 0.05 was considered statistically significant.

## References:

- [1] Development of a robust, low cost stem-loop real-time quantification PCR technique for miRNA expression analysis. *Mol Biol Rep.* 2013. PMID: 23307300.
- [2] Quantitative stem-loop RT-PCR for detection of microRNAs. *Methods Mol Biol.* 2011. PMID: 21533691.
- [3] MicroRNA detection in prostate tumors by quantitative real-time PCR (qPCR). *J Vis Exp.* 2012. PMID: 22643910.
- [4] Circulating microRNA-196a as a candidate diagnostic biomarker for chronic hepatitis C. *Mol Med Rep.* 2015. PMID: 25738504.
- [5] Serum microRNA-21 as marker for necroinflammation in hepatitis C patients with and without hepatocellular carcinoma. *PLoS One.* 2011. PMID: 22066022.
- [6] Diagnostic microRNA markers to screen for sporadic human colon cancer in stool: I. Proof of principle. *Cancer Genomics Proteomics.* 2013. PMID: 23741026.
- [7] Stem-loop RT-qPCR for microRNA expression profiling. *Methods Mol Biol.* 2012. PMID: 22144190.
- [8] Novel real-time PCR assay of microRNAs using S-Poly(T), a specific oligo(dT) reverse transcription primer with excellent sensitivity and specificity. *PLoS One.* 2012. PMID: 23152780.